



750 Garcia Ave, Pittsburg
CA 94565
+1.925.427.7800
HSIheadwalls.com

Bed Shoe Bed Locator

The floor mounted Bed Shoe Bed Locator will allow docking of one caster at the head end of the bed to assure proper location of the bed in relation to the Patient Service Column. The BedStop will allow the bed to pivot in an arc from the docked caster without losing its relationship to the column.

Description	Part Number
Bed Shoe	090800



Features/Specifications:

- Easy to install
- Durable
- Compact, space-saving solution